

NEW - SOMNOscreen™ plus PSG Tele BT:

Built-in Data Transmitter & Signal Control in the Patient's Room

Benefits at a glance

- Improved comfort with radio transmitter built-in to main unit
- Data Transfer from the Receiver to your network via LAN
- Optimized energy consumption, recording duration 30% longer
- Easy, fast and secure control of signals on a tablet PC
- Displays connect status with LED's
- All existing SOMNOscreen™ sensors are compatible
- No interference with existing measurement stations
- Lower potential for interference thanks to spatially adapted operating distance
- Modern color design



Art.Nr.: SBT202

Available as a set with Bluetooth IR Camera & built-in radio receiver

Contains: SOMNOscreen™ plus PSG Tele BT and Camera with built-in radio receiver

- Space-saving, dimensions: 135 x 105 x 60 mm
- Minimal installation efforts, only 1 LAN and 1 Power connection
- Minimum sensitivity to interference
- Resolution: 640 x 480 pixels with integrated IR-Illumination
- Allows Communication with the patient through built-in speaker and a microphone




Set with SOMNOscreen™ plus PSG Tele BT and Bluetooth-Camera VID100

Art.Nr.: SBTV202

Signal check easy, fast and secure on tablet PC

- Instant Signal Check on the patient's bedside
- Simply select the SOMNOscreen™ you wish to view from the tablet PC menu and it will connect to instantly display the patient's data



Room 1
Patient's Data 1

Room 2
Patient's Data 2

Room 3
Patient's Data 3