



Product Catalogue








SOMNOtouch™ NIBP











SOMNOtouch™ NIBP: Sets

Item no.:	Picture:	Description:
<p>TON100</p>		<p>SOMNOtouch NIBP Continuous non-invasive BP Recording</p> <p>Max. Recording time: 24 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse, Plethysmography, 3-channel-ECG, syst./diast. Blood pressure Patient marker</p> <p>Items included: Basic device SOMNOtouch NIBP, SpO2 Finger clip - either soft silicone or FlexiWrap (Adult), 3-channel-ECG, 1 pack ECG adhesive electrodes (50 pcs.), Nuprep ECG & EEG Abrasive Skin Prepping Gel, Watch strap size M, Dockingstation (incl. power supply and USB cable), Bag, User manual, Software DOMINO light (2 licences)</p>
<p>TON100A</p>		<p>Add-on device: SOMNOtouch NIBP Continuous non-invasive BP Recording</p> <p>Max. Recording time: 24 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse, Plethysmography, 3-channel-ECG, syst./diast. Blood pressure</p> <p>Items included: Basic device SOMNOtouch NIBP, SpO2 Finger clip - either soft silicone or FlexiWrap (Adult), 3-channel-ECG, 1 pack ECG adhesive electrodes (50 pcs.), Nuprep ECG & EEG Abrasive Skin Prepping Gel, Watch strap size M, Bag</p>





SOMNOtouch™ NIBP: Extended Analysis

Item no.:	Picture:	Description:
TOS031		Option for automatic Blood pressure calibration Items included: Pressure sensor for automatic Blood pressure calibration, Device for blood pressure measurements
TOS051		Simplified Sleep Stage Evaluation (ambulatory) Items included: EEG / EOG Combi sensor (3 channels) with integrated light sensor, Activation of EEG option (based on serial number), 1 pack of disposable electrodes
TOS100		Flow pressure option
TOS120		Analogue Optocoupler Connecting external devices to the SOMNOtouch Requires technical data sheet of the external device! (An additional charge for the adaption depending on the needed equipment will be added)
KOMP112B		Option Signal Check via Tablet PC (Bluetooth) including SOMNOmedics Android Application 1.6 GHz dual-core processor, 1 GB of memory, 16 GB internal memory, expandable via microSD (64GB) 10.1 inches (25.7 cm) multi-touch display (1280x800) WLAN, UMTS, Bluetooth 4.0, SIM card slot, Android 4.2.2, 2 years warranty (to ensure the functionality of the live data transmission, a SOMNOscreen BT Series or a SOMNOtouch is needed)
KOM112B		Option signal check via smartphone including SOMNOmedics Android application Signal check via e-mail for ambulatory polygraphy and polysomnography
SOW105		DOMINO light Zusatz-Softwarelizenz
SW010		Schiller Darwin OEM, ECG-Arrhythmia Analysis for DOMINO (light)






SOMNOtouch™ NIBP: Sensors











Item no:	Picture:	Description:
TOS010		3-channel ECG sensor for SOMNOtouch RESP and NIBP
TOS030		Pressure sensor for automatic Blood pressure calibration
TOS061		Soft Silicone Finger Sensor for SOMNOtouch NIBP 27,5 cm, for adults
TOS071		SoftWrap Finger Sensor for SOMNOtouch NIBP 27,5 cm, for adults
TOS050		EEG / EOG Combi sensor (3 channels) with integrated light sensor, for SOMNOtouch RESP
TOS902		Docking station for SOMNOtouch incl. power supply and USB cable
TOS903		3-port Docking station for SOMNOtouch incl. power supply and USB cable
SOW080		SOMNOmedics™ Bag for SOMNOtouch™

SOMNOtouch™ NIBP: Wristbands

Item no:	Picture:	Description:
TOS852		Watch Strap for SOMNOtouch NIBP circumference: ca. 18 - 23 cm, incl. 2 spring pins Silicone
TOS860		Watch Strap for SOMNOtouch NIBP Size S, circumference: ca. 13 - 20 cm fleece with velcro fastener, black
TOS861		Watch Strap for SOMNOtouch NIBP Size M, circumference: ca. 19 - 23 cm fleece with velcro fastener, black
TOS862		Watch Strap for SOMNOtouch NIBP Size L, circumference: ca. 23 - 30 cm fleece with velcro fastener, black




SOMNOtouch™ RESP: Electrodes and consumables

Item no:	Picture:	Description:
SEN002		Single use electrode N (recommended for ECG) Single use electrode with decentral snap fastener, 25 pcs./Pack
SEN004		Adhesive Electrode (recommended for ECG) Skin-friendly Santara-microporous fleece textile, 30 pcs./Pack
SEN010C		Disposable Electrodes (recommended for EOG, EEG) Silver-Silver chloride, self adhesive, Snap fastener, Diameter 2.5 cm, 30 pcs./Pack
SEN011		Disposable Electrodes Silver-Silver chloride, self adhesive, Snap fastener, Diameter 3 cm, 30 pcs./Pack
SEN012		Disposable Electrodes (recommended for PLM) Silver-Silver chloride, self adhesive, Snap fastener, 2.8 x 4.5 cm, 30 pcs./Pack

Item no:	Picture:	Description:
TOS075		Wrap tapes for SoftWrap Fingerclip 12 pcs./Box Single to use, highly flexible and ultra soft wrap tape
SEN039		Cleaning Pads, 10 pcs.
SEN115		Nasal Cannula HC (High Comfort) Adult, Nasal Pressure transducer with LuerLock connector male, Tube length 60 cm, 10 pcs./Pack
SEN116		Nasal Cannula Pediatric Nasal Pressure transducer with LuerLock connector male, Tube length 210 cm, 25 pcs./Pack
SEN117		Nasal Cannula Neonatal Nasal Pressure transducer with LuerLock connector male, Tube length 210 cm, 10 pcs./Pack
SEN118		Nasal/Oral Cannula for Pressure Sensor Adult Nasal and Oral Pressure transducer with LuerLock connector male, Tube length 60 cm, 25 pcs./Pack
SEN121		Nasal Cannula (divided) Adult Divided Nasal Pressure monitoring & Gas sampling cannula, Tube length 210 cm, 25 pcs./Pack
SEN903		Nasal Cannula (CO2) Adult for Capnometry, Tube length 210 cm, 10 pcs./Pack
SEN921		Silicone Tube for CPAP connection with LuerLock
SEN924		CPAP-Adapter



PCs

Item no:	Picture:	Description:
KOMP34W7C		PC Recording and analysis station min. Intel Core i5 3.2 GHz, 2x4 GB of memory, 1000 GB hard drive (7200 rev / min), Graphics Card: Intel HD Grafik4600 (1920x1200), 19 in 1 media card reader, 16 x DVD / - RW, Gigabit LAN, keyboard, mouse, Operating System Microsoft Windows 7 Ultimate 64 Bit OEM, Microsoft Office 2013 (Word, Excel, PowerPoint, Outlook) OEM, 3 years on-site service by Dell
KOMP83B		27 Inch LED-Flat Screen, resolution: 4k/3.840 x 2.160, incl. Soundbar
OPT011n		Duplex Laser Printer, 23 ppm color / 23 SW, 600x600dpi, Network port, USB 2.0, 256 MB, Paper cassette holds up to 250 sheets, max. monthly printing capacity of 40,000 pages