

Price list





SOMNOtouch™ RESP





Enduser International

Valid from 1st of September 2014















SOMNOtouch™ RESP: Sets

Item no.:	Picture:	Description:
TOR100		<p>SOMNOtouch™ RESP Respiratory Screening</p> <p>Max. Recording time: 19 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse rate, Plethysmography, 2x Effort, Flow+Snoring (Nasal pressure), CPAP Pressure, Patient marker</p> <p>Items included: Basic device, Belts(Abdomen/Thorax L), Effort sensor, SpO2 Finger clip Soft Silicone (Adult), 1 pack Nasal cannulas (10 pcs.), CPAP adapter, Dockingstation (incl. power supply and USB cable), SOMNOmedics™ Bag, User manual, Software DOMINO light (2 licences)</p>
TOR100A		<p>Add-on device: SOMNOtouch™ RESP Respiratory Screening</p> <p>Max. Recording time: 19 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse, Plethysmography, 2x Effort, Flow+Snoring (Nasal pressure), CPAP Pressure, Patient marker</p> <p>Items included: Basic device, Belts (Abdomen/Thorax L), SpO2 Finger clip Soft Silicone (Adult), 1 pack Nasal cannulas (10 pcs.), CPAP adapter, SOMNOmedics™ Bag</p>
TOR110		<p>SOMNOtouch™ RESP (RIP) Respiratory Screening</p> <p>Max. Recording time: 18 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse, Plethysmography, 2x RIP Effort, Flow, CPAP Pressure, Snoring (Nasal pressure), Patient marker</p> <p>Items included: Basic device, 2 Inductive effort belts Abdomen/Thorax L + Body mount, inductive effort sensor (RIP), SpO2 Finger clip Soft Silicone (Adult), 1 pack Nasal cannulas (10 pcs.), CPAP adapter, Dockingstation (incl. power supply and USB cable), SOMNOmedics™ Bag, User manual, Software DOMINO light (2 licences)</p>
TOR110A		<p>Add-on device: SOMNOtouch™ RESP (RIP) Respiratory Screening</p> <p>Max. Recording time: 18 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse rate, Plethysmography, 2x RIP Effort, Flow, CPAP Pressure, Snoring (Nasal pressure), Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, 2 Inductive effort belts Abdomen/Thorax size L + Body mount, Inductive effort sensor (RIP), SpO2 Finger clip Soft Silicone (Adult), 1 pack Nasal cannulas (10 pcs.), CPAP adapter, SOMNOmedics™ Bag</p>



Item no.:	Picture:	Description:
TOR105		<p>SOMNOtouch™ RESP (pediatric) Respiratory Screening</p> <p>Max. Recording time: 19 h</p> <p>11 channels:Body position, Movement, SpO2, Pulse rate, Plethysmography,2x Effort, Flow+Snoring (Nasal pressure), CPAP Pressure, Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, Belt set Abdomen/Thorax size XS, Effort sensor, SpO2 Finger clip Soft Silicone (pediatric), 1 pack Nasal cannulas ped. (10 pcs.), CPAP adapter, Dockingstation (incl. power supply and USB cable), SOMNOmedics™ Bag, manual, Software DOMINO light (2 licences)</p>
TOR105A		<p>Add-on device : SOMNOtouch™ RESP (pediatric) Respiratory Screening</p> <p>Max. Recording time: 19 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse rate, Plethysmography, 2x Effort, Flow+Snoring (Nasal pressure), CPAP Pressure, Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, Belt set Abdomen/Thorax size XS, Effort sensor, SpO2 Finger clip Soft Silicone (pediatric), 1 pack Nasal cannulas ped. (10 pcs.), CPAP adapter, SOMNOmedics™ Bag,</p>
TOR115		<p>SOMNOtouch RESP (RIP, pediatric) Respiratory Screening</p> <p>Max. Recording time: 18 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse, Plethysmography, 2x RIP Effort, Flow, CPAP Pressure, Snoring (Nasal pressure),Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, 2 Ind. effort belts Abdomen/Thorax size S + Body mount, Ind. effort sensor (RIP), SpO2 Finger clip Soft Silicone (pediatric), 1 pack Nasal cannulas ped. (10 pcs.), CPAP adapter, Dockingstation (incl. power supply and USB cable), SOMNOmedics™ Bag, manual, Software DOMINO light (2 licences)</p>
TOR115A		<p>Add-on device: SOMNOtouch RESP (RIP, pediatric) Respiratory Screening</p> <p>Max. Recording time: 18 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse, Pleth,2x RIP Effort, Flow, CPAP Pressure, Snoring (Nasal pressure),Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, 2 Inductive effort belts Abdomen/Thorax size S + Body mount, Inductive effort sensor (RIP), SpO2 Finger clip Soft Silicone (pediatric), 1 pack Nasal cannulas pediatric (10 pcs.), CPAP adapter, SOMNOmedics™ Bag</p>










SOMNOtouch™ RESP: Options for extended analysis







Item no.:	Picture:	Description:
TOS041		<p>1-channel EEG option for SOMNOtouch RESP</p> <p>Items included: 1-channel EEG sensor, Activation of EEG option (based on serial number), 1 pack of disposable electrodes</p>
TOS051		<p>Simplified Sleep Stage Evaluation (ambulatory) for SOMNOtouch RESP</p> <p>Items included: EEG / EOG Combi sensor (3 channels) with integrated light sensor, Activation of EEG option (based on serial number), 1 pack of disposable electrodes</p>
TOS095		<p>AASM Headbox for SOMNOtouch RESP</p> <p>Items included: AASM Headbox (4 EEG, 2 EOG, 1 EMG/ECG, Impedance check), Activation of AASM Sleep Option (based on serial number), 10 GRASS gold plated EEG cup electrodes, EC2 Electrode cream (self-adhesive), NUPREP Abrasive skin prepping gel</p>
TOS097		<p>R&K Headbox for SOMNOtouch RESP</p> <p>Items included: R&K Headbox (4 EEG, 2 EOG, 1 EMG/ECG, Impedance check), Activation of R&K Sleep Option (based on serial number), 10 GRASS gold plated EEG cup electrodes, EC2 Electrode cream (self-adhesive), NUPREP Abrasive skin prepping gel</p>
TOS091		<p>6-channel EEG option for SOMNOtouch RESP</p> <p>Items included: 6-channel EEG module (6 EEG, 1 EMG/ECG, Impedance control), Activation of EEG option (based on serial number), 10 GRASS gold plated EEG cup electrodes (10 mm with hole, Cable length: 76 cm), EC2 Electrode cream (self-adhesive), NUPREP Abrasive skin prepping gel</p>
VID98A		<p>LAN Infrared Video Camera Home Sleep</p> <p>H.264 Video Stream, Integrated IR illumination, Microphone, Resolution 800x600 pixels, Power supply, USB stick 8 GB for data storage, tripod, Software for synchronized recording and archiving</p>
TOS301		<p>Blood Pressure Option for SOMNOtouch RESP</p> <p>Items included: 1-channel ECG sensor for SOMNOtouch RESP, Activation of BP option (based on serial number), 1 Pack of Adhesive Electrodes</p>

Item no.:	Picture:	Description:
KOMP112C		Option Signal Check via Tablet PC (Bluetooth) including SOMNOmedics Android Application
KOMP111		Option Signal Check via smartphone including SOMNOmedics Android application Signal check via e-mail for ambulatory polygraphy and polysomnography
SSC091		DOMINO light Software option: extended Data Base for measurement results
SW010	 SCHILLER The Art of Diagnostics	Schiller Darwin OEM, ECG-Arrhythmia-Analysis for DOMINO light
SOW105		DOMINO light Additional Software License




SOMNOtouch™ RESP: Sensors/Accessories


Item no.:	Picture:	Description:
TOS011		1-channel ECG sensor for SOMNOtouch RESP
TOS012		1-channel ECG sensor pediatric for SOMNOtouch RESP

Item no.:	Picture:	Description:
TOS020		Effort Sensor für SOMNOtouch RESP
TOS060		Soft Silicone Finger Sensor for SOMNOtouch RESP 1.3 m, for adults
TOS062		Soft Silicone Finger Sensor for SOMNOtouch RESP 1.0 m, pediatric
TOS070		SoftWrap Finger Sensor for SOMNOtouch RESP 1.3 m, for adults and pediatric use
TOS080		Inductive effort sensor (RIP) for SOMNOtouch RESP 2 channels (Thorax/Abdomen)
TOS081		Set for inductive effort measurement (RIP) for SOMNOtouch RESP Items included: Inductive effort sensor 2 channels (Thorax/Abdomen), Inductive effort belts (2 pieces, size L), Body Mount SOMNOtouch for inductive effort belts
TOS083		Set for inductive effort measurement (RIP) with disposable belts for SOMNOtouch RESP Items included: Inductive effort sensor 2 channels (Thorax/Abdomen), Lock for inductive disposable effort belts (2x), Inductive disposable effort belts (20 m roll), Body Mount SOMNOtouch for inductive disposable effort belts
TOS110		2-channel PLM sensor for SOMNOtouch RESP
TOS120		Analogue Optocoupler Connecting external devices to the SOMNOtouch (An additional charge for the adaption depending on the needed equipment will be added)




Item no.:	Picture:	Description:
TOS140		Effort + 1PLM Combi Sensor for SOMNOtouch RESP
TOS151		Pneumotachograph Set for SOMNOtouch RESP Items included: a tube for the connection of CPAP/BiPAP devices, analyses the quantitative breathing and leakage flow
TOS170		Effort + 2PLM Combi Sensor for SOMNOtouch RESP
TOS902		Docking station for SOMNOtouch incl. power supply and USB cable
TOS903		3-port Docking station for SOMNOtouch incl. power supply and USB cable
SOW080		SOMNOmedics™ Bag for SOMNOwatch plus and SOMNOtouch

SOMNOtouch™ RESP: Belts






Item no.:	Picture:	Description:
TOS799		Effort belt for effort sensor SOMNOtouch RESP, size XS, circumference: ca. 50 - 75 cm, fleece with quick release and velcro fastener
TOS800		Effort belt for effort sensor SOMNOtouch RESP size S, circumference: ca. 65 - 88 cm, fleece with quick release and velcro fastener
TOS801		Effort belt for effort sensor SOMNOtouch RESP size M, circumference: ca. 74 - 109 cm, fleece with quick release and velcro fastener

Item no.:	Picture:	Description:
TOS802		Effort belt for effort sensor SOMNOtouch RESP size L, circumference: ca. 86 - 134 cm, fleece with quick release and velcro fastener
TOS803		Effort belt for effort sensor SOMNOtouch RESP size XXL, circumference: ca. 99 - 159 cm, fleece with quick release and velcro fastener
TOS804		Effort belt for SOMNOtouch RESP, size XS, Circumference: ca. 50 - 75 cm fleece with quick release and velcro fastener
TOS805		Effort belt for SOMNOtouch RESP size S, Circumference: ca. 65 - 88 cm, fleece with quick release and velcro fastener
TOS806		Effort belt for SOMNOtouch RESP size M, Circumference: ca. 74 - 109 cm fleece with quick release and velcro fastener
TOS807		Effort belt for SOMNOtouch RESP size L, Circumference: ca. 86 - 134 cm, fleece with quick release and velcro fastener
TOS808		Effort belt for SOMNOtouch RESP size XXL, Circumference: ca. 99 - 159 cm, fleece with quick release and velcro fastener
TOS810		Body Mount SOMNOtouch RESP, for inductive effort belts, Width 38 mm, Length 15 cm
TOS811		Application Mount SOMNOtouch RESP Requires effort belt TOS799 - TOS803
TOS812		Body Mount SOMNOtouch RESP, for inductive disposable effort belts, Width 38 mm, Length 18 cm
SEN660		Inductive Effort Belts - Size S Pediatric, (2 pcs./Pack), Circumference: ca. 40 - 70 cm
SEN661		Inductive Effort Belts - Size M Pediatric and Adult, (2 pcs./Pack), Circumference: ca. 60 - 90 cm
SEN662		Inductive Effort Belts - Size L Adult, (2 pcs./Pack), Circumference: ca. 80 - 120 cm
SEN663		Inductive Effort Belts - Size XXL Adult, (2 pcs./Pack), Circumference: ca. 100 - 150 cm

SOMNOtouch™ RESP: Creams








Item no.:	Picture:	Description:
SEN030		Nuprep Abrasive skin prepping gel, 114 g per tube, 3 pcs./Pack
SEN032		Ten20 EEG Conductive Adhesive Paste, 114 g tube, 3 pcs./Pack
SEN034		EC2 Adhesive electrode cream, water-based, 100 g per tube, 10 pcs./Pack

SOMNOtouch™ RESP: EEG-Electrodes

Item no.:	Picture:	Description:
SEN018		Connecting cables for disposable electrodes Snap fastener electrode cables (Ø 1.5 mm Safety pins), Cable length 80 cm, assorted colours, 10 pcs./Pack
F-E5GH-24		GRASS Gold Cup Electrodes 10 mm with hole 2 mm, assorted colours, Cable length 61 cm, 10 pcs./Box
F-E5GH-30		GRASS Gold Cup Electrodes 10 mm with hole 2 mm, assorted colours, Cable length 76 cm, 10 pcs./Box
F-E5GH-48		GRASS Gold Cup Electrodes 10 mm with hole 2 mm, assorted colours, Cable length 122 cm, 10 pcs./Box
F-E5SH-48		GRASS Silver Cup Electrodes 10 mm with hole 2 mm, assorted colours, Cable length 122 cm, 10 pcs./Box

SOMNOtouch™ RESP: Electrodes and consumables

Item no.:	Picture:	Description:
SEN002		Single use electrode N (recommended for ECG) Single use electrode with decentral snap fastener, 25 pcs./Pack
SEN004		Adhesive Electrode (recommended for ECG) Skin-friendly Santara-microporous fleece textile, 30 pcs./Pack
SEN010C		Disposable Electrodes (recommended for EOG, EEG) Silver-Silver chloride, self-adhesive, Snap fastener, Diameter 2.5 cm, 30 pcs./Pack
SEN011		Disposable Electrodes Silver-Silver chloride, self adhesive, Snap fastener, Diameter 3 cm, 30 pcs./Pack
SEN012		Disposable Electrodes (recommended for PLM and ECG) Silver-Silver chloride, self-adhesive, Snap fastener, 2.8 x 4.5 cm, 30 pcs./Pack
SEN020		Disposable Electrode with connecting cable Clear gel and Safety plug (colour coded) 1.5 mm, Cable length 60 cm, Diameter 2.2 cm, 15 pcs./Pack
SEN020A		Disposable Electrode with connecting cable Clear gel and Safety plug (colour coded) 1.5 mm, Cable length 120 cm, Diameter 2.2 cm, 15 pcs./Pack
SEN020B		Disposable Electrode with connecting cable Clear gel and Safety plug (colour coded) 1.5 mm, Cable length 160 cm, Diameter 2.2 cm, 15 pcs./Pack
TOS075		Wrap tapes for SoftWrap Fingerclip 12 pcs./Box, Single to use, highly flexible and ultra-soft wrap tape

Item no.:	Picture:	Description:
SEN115		Nasal Cannula HC (High Comfort) Adult, Nasal Pressure transducer with LuerLock connector male, Tube length 60 cm, 10 pcs./Pack
SEN115A		Nasal Cannula, Adult Nasal Pressure transducer with male LuerLock connector, Tube length 60 cm, 100 pcs./Pack
SEN116B		Nasal Cannula Pediatric Nasal Pressure transducer with LuerLock connector male, Tube length 40 cm, 10 pcs./Pack
SEN118		Nasal/Oral Cannula for Pressure Sensor Adult, Nasal and Oral Pressure transducer with LuerLock connector male, Tube length 60 cm, 25 pcs./Pack
SEN121		Nasal Cannula (divided) Adult divided Nasal Pressure monitoring & Gas sampling cannula, Tube length 210 cm, 25 pcs./Pack
SEN921		Silicone Tube for CPAP connection with LuerLock
SEN924		CPAP adapter

PCs/Notebooks are available in different configurations on request***